

DyLight 488, Goat Anti-Mouse IgG

Description

Product type	Secondary Antibody
Code	AS0806
Host	Goat
Isotype	N/A
Size	100ul, 1ml
Immunogen	N/A
Mol wt	N/A
Species reactivity	Mouse
Clonality	Polyclonal
Recommended application	IF, ICC, FCM
Concentration	1mg/mL
Full name	DyLight 488, Goat Anti-Mouse IgG
Synonyms	N/A

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Secondary antibodies are available conjugated to enzyme, biotin or fluorophore for use in a variety of antibody-based applications including Western Blot, ImmunoHistoChemistry, ImmunoFluorescence, Flow Cytometry and ELISA. We offer high quality secondary antibodies from goat, rabbit and donkey sources for your each application. Serum adsorbed secondary antibodies are also available and are recommended for use with immunoglobulin-rich samples.

Recommended Dilution

Not yet tested in other applications.

Images

No images.

Storage

-20°C for one year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com