

A-FABP Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP00166
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human A-FABP
Mol wt	14520
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	A-FABP
Synonyms	FABP4; 422/aP2; A-FABP; adipocyte; ALBP; aP2; Fatty acid binding protein 4 adipocyte; HEL S 104; Fatty acid-binding protein; LBP1

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

FABP4 encodes the fatty acid binding protein found in adipocytes. Fatty acid binding proteins are a family of small| highly conserved| cytoplasmic proteins that bind long-chain fatty acids and other hydrophobic ligands. It is thought that FABPs roles include fatty acid uptake| transport| and metabolism.

Recommended Dilution

WB: 1: 1000 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com