

## ADAMTS-7 Polyclonal Antibody

### Description

<b>Product type</b>	Primary Antibody
<b>Code</b>	BT-AP00250
<b>Host</b>	Rabbit
<b>Isotype</b>	IgG
<b>Size</b>	20ul, 50ul, 100ul
<b>Immunogen</b>	Synthesized peptide derived from ADAMTS-7 . at AA range: 150-230
<b>Mol wt</b>	109695
<b>Species reactivity</b>	Human
<b>Clonality</b>	Polyclonal
<b>Recommended application</b>	IHC-p, ELISA
<b>Concentration</b>	1 mg/ml
<b>Full name</b>	ADAMTS-7 Antibody
<b>Synonyms</b>	ADAMTS7; A disintegrin and metalloproteinase with thrombospondin motifs 7; ADAM-TS 7; ADAM-TS7; ADAMTS-7; COMPase

**This product is for research use only, not for use in human, therapeutic or diagnostic procedure.**

### Background

The ADAM metalloproteinase with thrombospondin type 1 motif 7 encoded by ADAMTS7 is a member of the ADAMTS (a disintegrin and metalloproteinase with thrombospondin motifs) family. Members of this family share several distinct protein modules, including a propeptide region, a metalloproteinase domain, a disintegrin-like domain, and a thrombospondin type 1 (TS) motif. Individual members of this family differ in the number of C-terminal TS motifs, and some have unique C-terminal domains. The encoded preproprotein is proteolytically processed to generate the mature enzyme. This enzyme contains two C-terminal TS motifs and may regulate vascular smooth muscle cell (VSMC) migration. Mutations in this gene may be associated with susceptibility to coronary artery disease.

### Recommended Dilution

IHC: 1: 100 - 1: 300

ELISA: 1: 40000

Not yet tested in other applications.

### Images

No images.

### Storage

-20°C for one year