

Cdk9 (Acetyl Lys44) Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP01333
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human Cdk9 (Acetyl Lys44)
Mol wt	40920
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	Cdk9
Synonyms	Cdk9 ;Acetyl Lys44; Cyclin-dependent kinase 9; EC 2.7.11.22; EC 2.7.11.23; C-2K; Cell division cycle 2-like protein kinase 4; Cell division protein kinase 9; Serine/threonine-protein kinase PITALRE; Tat-associated kinase complex catalytic subunit

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

The protein encoded by this gene is a member of the cyclin-dependent protein kinase (CDK) family. CDK family members are highly similar to the gene products of *S. cerevisiae* cdc28 and *S. pombe* cdc2 and known as important cell cycle regulators. This kinase was found to be a component of the multiprotein complex TAK/P-TEFb which is an elongation factor for RNA polymerase II-directed transcription and functions by phosphorylating the C-terminal domain of the largest subunit of RNA polymerase II. This protein forms a complex with and is regulated by its regulatory subunit cyclin T or cyclin K. HIV-1 Tat protein was found to interact with this protein and cyclin T which suggested a possible involvement of this protein in AIDS.

Recommended Dilution

WB: 1: 1000 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year