

CUTC Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP02033
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human CUTC
Mol wt	30030
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	Copper homeostasis protein cutC homolog
Synonyms	Copper homeostasis protein cutC homolog

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Recommended Dilution

WB: 1: 1000 - 1: 2000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com