

DAP-1 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP02496
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	The antiserum was produced against synthesized peptide derived from human DAP. AA range:11-60
Mol wt	11165
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	IHC-p, ELISA
Concentration	1 mg/ml
Full name	DAP-1 Antibody
Synonyms	DAP; DAP1; Death-associated protein 1; DAP-1

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

DAP encodes a basic, proline-rich, 15-kD protein. Death-associated protein 1 acts as a positive mediator of programmed cell death that is induced by interferon-gamma. Alternatively spliced transcript variants encoding distinct isoforms have been found for this gene.


Recommended Dilution

IHC: 1: 100 - 1: 300

ELISA: 1: 40000

Not yet tested in other applications.

Images



Immunohistochemistry analysis of paraffin-embedded human colon carcinoma tissue, using DAP Antibody. The picture on the right is blocked with the synthesized peptide.

Storage

-20°C for one year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com