

## FA18A Rabbit Polyclonal Antibody

### Description

Product type	Primary Antibody
Code	BT-AP02568
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	Synthesized peptide derived from human FA18A
Mol wt	23430
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	FA18A
Synonyms	FA18A

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

### Background


This gene encodes a membrane protein associated with the Golgi apparatus| which plays a crucial role in intracellular vesicular transport. The encoded protein is likely associated with the late (trans) Golgi compartments| which are involved in the delivery of secretory and membrane proteins to the endosome| lysosome or the plasma membrane. Alternative splicing results in multiple transcript variants.

### Recommended Dilution

WB: 1: 500 - 1: 2000

Not yet tested in other applications.

### Images



### Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: [save@bt-laboratory.com](mailto:save@bt-laboratory.com) | [www.bt-laboratory.com](http://www.bt-laboratory.com)