

FA8B Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP02611
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human FA8B (Cleaved-Ser760)
Mol wt	258610
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	FA8B
Synonyms	FA8B ;Cleaved-Ser760; Coagulation factor VIII; Antihemophilic factor; AHF; Procoagulant component; Factor VIIIa heavy chain; 200 kDa isoform; Factor VIIIa heavy chain; 92 kDa isoform; Factor VIII B chain; Factor VIIIa light chain;

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

disease:Defects in F8 are the cause of hemophilia A (HEMA)

Recommended Dilution

WB: 1: 1000 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com