

## FoxO1A (Acetyl Lys245) Rabbit Polyclonal Antibody

### Description

<b>Product type</b>	Primary Antibody
<b>Code</b>	BT-AP02718
<b>Host</b>	Rabbit
<b>Isotype</b>	IgG
<b>Size</b>	100ul, 50ul, 20ul
<b>Immunogen</b>	Synthesized peptide derived from human FoxO1A (Acetyl Lys245)
<b>Mol wt</b>	72050
<b>Species reactivity</b>	Human, Mouse, Rat
<b>Clonality</b>	Polyclonal
<b>Recommended application</b>	WB, ELISA
<b>Concentration</b>	1 mg/ml
<b>Full name</b>	FoxO1A
<b>Synonyms</b>	FoxO1A ;Acetyl Lys245; Forkhead box protein O1; Forkhead box protein O1A; Forkhead in rhabdomyosarcoma

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

### Background

disease:Chromosomal aberrations involving FOXO1 are a cause of rhabdomyosarcoma 2 (RMS2)


### Recommended Dilution

WB: 1: 1000 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

### Images



Western Blot analysis of HEK 293 cells, cell treated or untreated by TSA 400nM 24h. Primary Antibody was diluted at 1:1000. Secondary antibody was diluted at 1:10000

### Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: [save@bt-laboratory.com](mailto:save@bt-laboratory.com) | [www.bt-laboratory.com](http://www.bt-laboratory.com)