

G3ST1 Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP02761
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human G3ST1
Mol wt	46530
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	Galactosylceramide sulfotransferase
Synonyms	Galactosylceramide sulfotransferase ;GalCer sulfotransferase;EC 2.8.2.11;3'-phosphoadenosine-5'-phosphosulfate:GalCer sulfotransferase;3'-phosphoadenylylsulfate:galactosylceramide 3'-sulfotransf

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

N/A

Recommended Dilution

WB: 1: 1000 - 1: 2000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com