

IL-7 Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP03084
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human IL-7
Mol wt	19470
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	IL-7
Synonyms	IL-7; Interleukin-7; IL-7

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Hematopoietic growth factor capable of stimulating the proliferation of lymphoid progenitors. It is important for proliferation during certain stages of B-cell maturation. [online information: Interleukin-7 entry] [Belongs to the IL-7/IL-9 family.]

Recommended Dilution

WB: 1: 1000 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com