

MAFA Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP03527
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human MAFA
Mol wt	N/A
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	MAFA
Synonyms	MAFA; Transcription factor MafA; Pancreatic beta-cell-specific transcriptional activator; Transcription factor RIPE3b1; V-maf musculoaponeurotic fibrosarcoma oncogene homolog A

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background


MAFA is a transcription factor that binds RIPE3b| a conserved enhancer element that regulates pancreatic beta cell-specific expression of the insulin gene (INS; MIM 176730) (Olbrot et al. | 2002

Recommended Dilution

WB: 1: 500 - 1: 2000

Not yet tested in other applications.

Images



Western blot analysis of lysates from MDA-MB cells, primary antibody was diluted at 1:1000, 4°C overnight

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com