

NBPFG Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP03771
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human NBPFG
Mol wt	73700
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	NBPFG
Synonyms	NBPFG

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background


This gene is a member of the neuroblastoma breakpoint family (NBPF) which consists of dozens of recently duplicated genes primarily located in segmental duplications on human chromosome 1. This gene family has experienced its greatest expansion within the human lineage and has expanded to a lesser extent among primates in general. Members of this gene family are characterized by tandemly repeated copies of DUF1220 protein domains. Gene copy number variations in the human chromosomal region 1q21.1 where most DUF1220 domains are located have been implicated in a number of developmental and neurogenetic diseases such as microcephaly| macrocephaly| autism| schizophrenia| mental retardation| congenital heart disease| neuroblastoma| and congenital kidney and urinary tract anomalies. Altered expression of some gene family members is associated with several types of cancer. This gene family contains numerous pseudogenes.

Recommended Dilution

WB: 1: 500 - 1: 2000

Not yet tested in other applications.

Images



Western blot analysis of lysates from HeLa cells, primary antibody was diluted at 1:1000, 4°C overnight

Storage

-20°C for 1 year