

NPS3A Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP03860
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human NPS3A
Mol wt	27170
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	NPS3A
Synonyms	NPS3A

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background


NIPSNAP3A belongs to a family of proteins with putative roles in vesicular transport (Buechler et al., 2004)

Recommended Dilution

WB: 1: 500 - 1: 2000

Not yet tested in other applications.

Images



Western blot analysis of lysates from MCF-7 cells, primary antibody was diluted at 1:1000, 4°C overnight

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com