

PDE1C Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP04081
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human PDE1C
Mol wt	77990
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	Calcium/calmodulin-dependent 3'
Synonyms	Calcium/calmodulin-dependent 3'; 5'-cyclic nucleotide phosphodiesterase 1C ;Cam-PDE 1C;EC 3.1.4.17;hCam-3

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Recommended Dilution

WB: 1: 1000 - 1: 2000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhui Dist, Jiading, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com