

HLA-H Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP04135
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	The antiserum was produced against synthesized peptide derived from human HLAH. AA range:71-120
Mol wt	40850
Species reactivity	Human
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	HLA-H Antibody
Synonyms	HLA-H; HLAH; Putative HLA class I histocompatibility antigen; alpha chain H; HLA-12.4; HLA-AR; MHC class I antigen H

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

This major histocompatibility complex gene represents a transcribed pseudogene, possibly derived from HLA-A. HLA-H displays extensive variation.



Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 10000

Not yet tested in other applications.

Images



Storage

-20°C for one year

