

## PI51A Rabbit Polyclonal Antibody

### Description

Product type	Primary Antibody
Code	BT-AP04148
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human PI51A
Mol wt	N/A
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	PI51A
Synonyms	PI51A; Phosphatidylinositol 4-phosphate 5-kinase type-1 alpha; PIP5K1-alpha; PtdIns;4P-5-kinase 1 alpha; EC 2.7.1.68; 68 kDa type I phosphatidylinositol 4-phosphate 5-kinase alpha; Phosphatidylinositol 4-phosphate 5-kinase type I alpha; PIP5K1alpha

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

### Background


N/A

### Recommended Dilution

WB: 1: 500 - 1: 2000

Not yet tested in other applications.

### Images



Western blot analysis of lysates from AD293 cells, primary antibody was diluted at 1:1000, 4°C overnight

### Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: [save@bt-laboratory.com](mailto:save@bt-laboratory.com) | [www.bt-laboratory.com](http://www.bt-laboratory.com)