

PZR(Phospho-Tyr263) Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP04352
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	Synthesized phosho peptide around human PZR (Tyr263)
Mol wt	N/A
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	PZR
Synonyms	PZR ;Tyr263 ; Myelin protein zero-like protein 1; Protein zero-related

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background


N/A

Recommended Dilution

WB: 1: 1000 - 1: 2000

Not yet tested in other applications.

Images



Western blot analysis of lysates from SH-SY5Y cells, primary antibody was diluted at 1:1000, 4°C overnight

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhui Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com