

SIRPD Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP04754
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human SIRPD
Mol wt	21670
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	Signal-regulatory protein delta
Synonyms	Signal-regulatory protein delta ;SIRP-delta;Protein tyrosine phosphatase non-receptor type substrate 1-like 2

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Recommended Dilution

WB: 1: 1000 - 1: 2000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year