

ST65G Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP04918
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human ST65G
Mol wt	45540
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	STAGA complex 65 subunit gamma
Synonyms	STAGA complex 65 subunit gamma ;Adenocarcinoma antigen ART1;SPTF-associated factor 65 gamma;STAF65gamma;Suppressor of Ty 7-like

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Recommended Dilution

WB: 1: 1000 - 1: 2000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhui Dist, Jiading, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com