

TM55A Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP05172
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	Synthesized peptide derived from human TM55A
Mol wt	28270
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	TM55A
Synonyms	TM55A

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background


TMEM55A catalyzes the degradation of phosphatidylinositol 4|5-bisphosphate (PtdIns-4|5-P2) by removing the 4-phosphate (Ungewickell et al. 2005)

Recommended Dilution

WB: 1: 500 - 1: 2000

Not yet tested in other applications.

Images



Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com