

TTL3 Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP05313
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	Synthesized peptide derived from human TTL3
Mol wt	N/A
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	TTL3
Synonyms	TTL3; Tubulin monoglycylase TTL3; EC 6.3.2.-; HOTTTL; Tubulin--tyrosine ligase-like protein 3

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.


Background

Recommended Dilution

WB: 1: 500 - 1: 2000

Not yet tested in other applications.

Images



Western blot analysis of lysates from HpeG2 cells, primary antibody was diluted at 1:1000, 4°C overnight

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com