

UQCC Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP05413
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	Synthesized peptide derived from human UQCC
Mol wt	32890
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	UQCC
Synonyms	UQCC

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background


This gene encodes a transmembrane protein that is structurally similar to the mouse basic fibroblast growth factor repressed ZIC-binding protein. In mouse this protein may be involved in fibroblast growth factor regulated growth control. In humans polymorphisms in this gene are associated with variation in human height and osteoarthritis. Alternate splicing results in multiple transcript variants.

Recommended Dilution

WB: 1: 500 - 1: 2000

Not yet tested in other applications.

Images



Western blot analysis of lysates from 293T cells, primary antibody was diluted at 1:1000, 4°C overnight

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com