

WHAMM Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP05507
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human WHAMM
Mol wt	88990
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	WASP homolog-associated protein with actin membranes and microtubules
Synonyms	WASP homolog-associated protein with actin; membranes and microtubules ;WAS protein homology region 2 domain-containing protein 1;WH2 domain-containing protein 1

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Recommended Dilution

WB: 1: 1000 - 1: 2000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year