

Ndfip1 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP05824
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	Synthesized peptide derived from Ndfip1 . at AA range: 150-230
Mol wt	24899
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	Ndfip1 Antibody
Synonyms	NDFIP1; N4WBP5; PSEC0192; PSEC0223; NEDD4 family-interacting protein 1; Breast cancer-associated protein SGA-1M; NEDD4 WW domain-binding protein 5; Putative MAPK-activating protein PM13; Putative NF-k

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

The protein encoded by NDFIP1 (Nedd4 family interacting protein 1) belongs to a small group of evolutionarily conserved proteins with three transmembrane domains. It is a potential target for ubiquitination by the Nedd4 family of proteins. This protein is thought to be part of a family of integral Golgi membrane proteins.


Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 40000

Not yet tested in other applications.

Images



Western blot analysis of lysates from Jarkat cells, primary antibody was diluted at 1:1000, 4°over night

Storage

-20°C for one year