

B3GN1 Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP06250
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human B3GN1
Mol wt	45650
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	N-acetyllactosaminide beta-1 3-N-acetylglucosaminyltransferase
Synonyms	N-acetyllactosaminide beta-1; 3-N-acetylglucosaminyltransferase; EC 2.4.1.149; I-beta-1; 3-N-acetylglucosaminyltransferase; iGnT; Poly-N-acetyllactosamine extension enzyme; UDP-GlcNAc:betaGal beta-1

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Recommended Dilution

WB: 1: 1000 - 1: 2000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com