

BFP Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP06749
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	Recombinant Protein of BFP
Mol wt	N/A
Species reactivity	Species independent
Clonality	Polyclonal
Recommended application	WB

Concentration

Full name

Synonyms

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Blue Fluorescent Protein (BFP) emits a strong blue fluorescence. EBFP can be used as protein "tags" to study the sub-cellular localization of proteins and/or their translocation upon stimulation or as markers for transfection in transient and stable expression systems.

Recommended Dilution

WB: 1: 5000

Not yet tested in other applications.

Images



Western blot analysis of EBFP recombinant protein, diluted at 1:5000. Secondary antibody was diluted at 1:20000

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com