

BFP Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP06749
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	Recombinant Protein of BFP
Mol wt	N/A
Species reactivity	Species independent
Clonality	Polyclonal
Recommended application	WB
Concentration	N/A
Full name	N/A
Synonyms	N/A

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Blue Fluorescent Protein (BFP) emits a strong blue fluorescence. EBFP can be used as protein "tags" to study the sub-cellular localization of proteins and/or their translocation upon stimulation or as markers for transfection in transient and stable expression systems.

Recommended Dilution

WB: 1: 5000

Not yet tested in other applications.

Images



Western blot analysis of EBFP recombinant protein, diluted at 1:5000. Secondary antibody was diluted at 1:20000

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com