

PAGE-3 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP06882
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	Synthesized peptide derived from PAGE-3 . at AA range: 40-120
Mol wt	12480
Species reactivity	Human
Clonality	Polyclonal
Recommended application	IHC-p, ELISA
Concentration	1 mg/ml
Full name	PAGE-3 Antibody
Synonyms	PAGE3; GAGED1; Putative G antigen family D member 1; Prostate-associated gene 3 protein; PAGE-3

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

PAGE3 is a member of family of proteins that are expressed in a variety of tumors and in some fetal and reproductive tissues. Multiple alternatively spliced transcript variants have been observed.

Recommended Dilution

IHC: 1: 100 - 1: 300

ELISA: 1: 5000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for one year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com