

RBM16 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP07710
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	Synthesized peptide derived from RBM16 . at AA range: 550-630
Mol wt	140519
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	RBM16 Antibody
Synonyms	SCAF8; CCAP7; KIAA1116; RBM16; Protein SCAF8; CDC5L complex-associated protein 7; RNA-binding motif protein 16; SR-related and CTD-associated factor 8

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

May play a role in mRNA processing.


Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 10000

Not yet tested in other applications.

Images



Western Blot analysis of JK cells using RBM16 Polyclonal Antibody cells nucleus.

Storage

-20°C for one year