

RFPL4B Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP07778
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	The antiserum was produced against synthesized peptide derived from human RFPL4B. AA range:100-149
Mol wt	29850
Species reactivity	Human
Clonality	Polyclonal
Recommended application	IHC-p, ELISA
Concentration	1 mg/ml
Full name	RFPL4B Antibody
Synonyms	RFPL4B; RNF211; Ret finger protein-like 4B; RING finger protein 211

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

RFPL4B (Ret Finger Protein Like 4B) is a Protein Coding gene. An important paralog of this gene is RFPL4A.

Recommended Dilution

IHC: 1: 100 - 1: 300

ELISA: 1: 40000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for one year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com