

CNTN3 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP07917
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from part region of human protein
Mol wt	N/A
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	Contactin-3
Synonyms	Contactin-3 ;Brain-derived immunoglobulin superfamily protein 1;BIG-1;Plasmacytoma-associated neuronal glycoprotein

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Contactins mediate cell surface interactions during nervous system development. Has some neurite outgrowth-promoting activity. Belongs to the immunoglobulin superfamily. Contactin family. Contains 4 fibronectin type-III domains. Contains 6 Ig-like C2-type (immunoglobulin-like) domains. tissue specificity: In brain, it is expressed in frontal lobe, occipital lobe, cerebellum and amygdala.

Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com