

Slit3 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP08364
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	Synthesized peptide derived from the Internal region of human Slit3.
Mol wt	167685
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	IHC-p, ELISA
Concentration	1 mg/ml
Full name	Slit3 Antibody
Synonyms	SLIT3; KIAA0814; MEGF5; SLIL2; Slit homolog 3 protein; Slit-3; Multiple epidermal growth factor-like domains protein 5; Multiple EGF-like domains protein 5

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

The protein encoded by SLIT3 (slit guidance ligand 3) is secreted, likely interacting with roundabout homolog receptors to effect cell migration. Two transcript variants encoding different isoforms have been found for SLIT3.

Recommended Dilution

IHC: 1: 100 - 1: 300

ELISA: 1: 40000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for one year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com