

DUS15 Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP09363
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	Synthesized peptide derived from human DUS15
Mol wt	32450
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	DUS15
Synonyms	DUS15

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background


The protein encoded by this gene has both protein-tyrosine phosphatase activity and serine/threonine-specific phosphatase activity| and therefore is known as a dual specificity phosphatase. This protein may function in the differentiation of oligodendrocytes. Alternative splicing results in multiple transcript variants.

Recommended Dilution

WB: 1: 500 - 1: 2000

Not yet tested in other applications.

Images



Western blot analysis of lysates from U2OS cells, primary antibody was diluted at 1:1000, 4°C overnight

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com