

EST1A Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP09763
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human EST1A
Mol wt	156090
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	Telomerase-binding protein EST1A
Synonyms	Telomerase-binding protein EST1A ;EC 3.1.-.-;EST1-like protein A;Ever shorter telomeres 1A;Smg-6 homolog;Telomerase subunit EST1A;hSmg5/7a

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Recommended Dilution

WB: 1: 1000 - 1: 2000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com