

His-Tag Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP09900
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthetic peptide:HHHHHH conjugated to KLH.
Mol wt	N/A
Species reactivity	Species independent
Clonality	Polyclonal
Recommended application	WB
Concentration	
Full name	His Tag
Synonyms	His Tag

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

A polyhistidine-tag is an amino acid motif in proteins that consists of at least five histidine (His) residues, often at the N- or C-terminus of the protein. It is also known as hexa histidine-tag, 6xHis-tag, and by the trademarked name His-tag. Polyhistidine-tags are often used for affinity purification of polyhistidine-tagged recombinant proteins expressed in Escherichia coli and other prokaryotic expression systems.

Recommended Dilution

WB: 1: 5000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year