

GBG11 Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP10130
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human GBG11
Mol wt	8030
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	IHC-p, IF
Concentration	1 mg/ml
Full name	GBG11
Synonyms	GBG11

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

This gene is a member of the guanine nucleotide-binding protein (G protein) gamma family and encodes a lipid-anchored cell membrane protein. As a member of the heterotrimeric G protein complex this protein plays a role in this transmembrane signaling system. This protein is also subject to carboxyl-terminal processing. Decreased expression of this gene is associated with splenic marginal zone lymphomas.

Recommended Dilution

IHC-p: 1: 50 - 1: 200

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com