

JAK3(Phospho Tyr981) Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP10726
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human JAK3 (Phospho Tyr981)
Mol wt	123640
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	JAK3
Synonyms	JAK3 ;Phospho Tyr981; Tyrosine-protein kinase JAK3; EC 2.7.10.2; Janus kinase 3; JAK-3; Leukocyte janus kinase; L-JAK

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Catalytic activity:ATP + a

Recommended Dilution

WB: 1: 1000 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com