

LEAP2 Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP11034
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human LEAP2
Mol wt	8470
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	IHC-p, IF
Concentration	1 mg/ml
Full name	LEAP2
Synonyms	LEAP2

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

This gene encodes a cysteine-rich cationic antimicrobial peptide that is expressed predominantly in the liver. The mature peptide has activity against gram-positive bacteria and yeasts.

Recommended Dilution

IHC-p: 1: 50 - 1: 200

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com