

MD12L Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP11276
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human MD12L
Mol wt	235950
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	MD12L
Synonyms	MD12L

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background


The protein encoded by this gene is part of the Mediator complex| which is involved in transcriptional coactivation of nearly all RNA polymerase II-dependent genes. The Mediator complex links gene-specific transcriptional activators with the basal transcription machinery.

Recommended Dilution

WB: 1: 500 - 1: 2000

Not yet tested in other applications.

Images



Western blot analysis of lysates from HeLa cells, primary antibody was diluted at 1:1000, 4°C overnight

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com