

NAT10 Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP11495
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human NAT10
Mol wt	112750
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	NAT10
Synonyms	NAT10

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background


The protein encoded by this gene is an RNA cytidine acetyltransferase involved in histone acetylation| tRNA acetylation| the biosynthesis of 18S rRNA| and the enhancement of nuclear architecture and chromatin organization.

Recommended Dilution

WB: 1: 500 - 1: 2000

Not yet tested in other applications.

Images



Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com