

mStrawberry Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP11556
Host	Mouse
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	Recombinant Protein of mStrawberry
Mol wt	N/A
Species reactivity	Species independent
Clonality	Polyclonal
Recommended application	WB
Concentration	
Full name	
Synonyms	

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

mStrawberry is a bright, monomeric red fluorescent protein. mStrawberry excitation and emission maxima are 574 and 596 nm, respectively.

Recommended Dilution

WB: 1: 5000

Not yet tested in other applications.

Images



Western blot analysis of mStrawberry recombinant protein, diluted at 1:5000. Secondary antibody was diluted at 1:20000

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com