

Myc-Tag Polyclonal Antibody, HRP Conjugated

Description

Product type	Primary Antibody
Code	BT-AP11602
Host	Rabbit
Isotype	IgG
Size	100ul
Immunogen	N/A
Mol wt	N/A
Species reactivity	Species independent
Clonality	Polyclonal
Recommended application	WB
Concentration	1mg/ml
Full name	Myc-Tag
Synonyms	Myc-Tag; Myc

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

N/A

Recommended Dilution

WB: 1: 5000

ELISA: 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com