

MYO1B Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP11640
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human protein . at AA range: 800-880
Mol wt	N/A
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	Unconventional myosin-Ib
Synonyms	Unconventional myosin-Ib ;MYH-1c;Myosin I alpha;MMI-alpha;MMIa

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Motor protein that may participate in process critical to neuronal development and function such as cell migration, neurite outgrowth and vesicular transport.,Contains 1 myosin head-like domain.,Contains 6 IQ domains.,

Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com