

p53 (Mono Methyl Lys372) Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP11712
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human p53 (Mono Methyl Lys372)
Mol wt	43230
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	p53
Synonyms	p53 ;Mono Methyl Lys372; Cellular tumor antigen p53; Antigen NY-CO-13; Phosphoprotein p53; Tumor suppressor p53

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

cofactor: Binds 1 zinc ion per subunit. | disease: Defects in TP53 are a cause of choroid plexus papilloma

Recommended Dilution

WB: 1: 1000 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com