

RN216 Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP12319
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human RN216
Mol wt	95260
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	E3 ubiquitin-protein ligase RNF216
Synonyms	E3 ubiquitin-protein ligase RNF216 ;EC 6.3.2.-;RING finger protein 216;Triad domain-containing protein 3;Ubiquitin-conjugating enzyme 7-interacting protein 1;Zinc finger protein inhibiting NF-

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Recommended Dilution

WB: 1: 1000 - 1: 2000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhui Dist, Jiading, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com