

S17A5 Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP12372
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human S17A5
Mol wt	54450
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	Sialin
Synonyms	Sialin ;Membrane glycoprotein HP59;Sodium/sialic acid cotransporter;AST;Solute carrier family 17 member 5

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

N/A

Recommended Dilution

WB: 1: 1000 - 1: 2000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year