

## S2535 Rabbit Polyclonal Antibody

### Description

Product type	Primary Antibody
Code	BT-AP12381
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	Synthesized peptide derived from human S2535
Mol wt	33000
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	S2535
Synonyms	S2535

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

### Background


SLC25A35 belongs to the SLC25 family of mitochondrial carrier proteins (Haitina et al. 2006)

### Recommended Dilution

WB: 1: 500 - 1: 2000

Not yet tested in other applications.

### Images



Western blot analysis of lysates from SW480 cells, primary antibody was diluted at 1:1000, 4°C overnight

### Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: [save@bt-laboratory.com](mailto:save@bt-laboratory.com) | [www.bt-laboratory.com](http://www.bt-laboratory.com)