

SIRT1(Phospho-Ser549) Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP12661
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human SIRT1 (phospho-Ser549)
Mol wt	82170
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	SIRT1
Synonyms	SIRT1 ;phospho-Ser549; NAD-dependent protein deacetylase sirtuin-1; hSIRT1; EC 3.5.1.-; Regulatory protein SIR2 homolog 1; SIR2-like protein 1; hSIR2; SirtT1 75 kDa fragment; 75SirT1;

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

This gene encodes a member of the sirtuin family of proteins| homologs to the yeast Sir2 protein. Members of the sirtuin family are characterized by a sirtuin core domain and grouped into four classes. The functions of human sirtuins have not yet been determined; however| yeast sirtuin proteins are known to regulate epigenetic gene silencing and suppress recombination of rDNA. Studies suggest that the human sirtuins may function as intracellular regulatory proteins with mono-ADP-ribosyltransferase activity. The protein encoded by this gene is included in class I of the sirtuin family. Alternative splicing results in multiple transcript variants.

Recommended Dilution

WB: 1: 1000 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com