

PM34 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP13194
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	Synthesized peptide derived from human protein . at AA range: 140-220
Mol wt	N/A
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	Peroxisomal membrane protein PMP34
Synonyms	Peroxisomal membrane protein PMP34 ;34 kDa peroxisomal membrane protein;Solute carrier family 25 member 17

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

This gene encodes a peroxisomal membrane protein that belongs to the family of mitochondrial solute carriers. It is expressed in the liver, and is likely involved in transport. Alternative splicing results in multiple transcript variants.


Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images



Western blot analysis of lysates from K562 cells, primary antibody was diluted at 1:1000, 4°C overnight

Storage

-20°C for 1 year